


Implementation Guide



mobilize
rescue systems

Equip. Empower. **Respond.**

Rescue Systems User Guide



Welcome to the Mobilize Community!

You are just a few steps away from empowering your entire organization to respond to medical emergencies.

1. How to install Public Access Wall Boxes and Alarmed Wall Cabinets

Mount Public Access Wall Boxes and Alarmed Wall Cabinets in high-traffic areas with clear visibility and easy access

Mounting and alarming instructions can be found in the manual links on the left of this page

Mobilize Recommends:

- Mount securely to wall
- Write licensing code on the Public Access Wall Boxes
- Set and test alarm wearing ear protection
- Read [ADA GUIDELINES](#)



Public Access Wall Boxes

[MODEL 50](#)
[MODEL 100](#)
[HOW TO ALARM](#)



Comprehensive & Mobile

[ALARMED WALL CABINET](#)

2. How to register and activate the Mobilize Rescue App in your units

COMPREHENSIVE & MOBILE UNITS:

ENSURE UNITS ARE FULLY CHARGED

- Connect to WiFi for the initial activation of each device
 - WiFi connection is in the settings menu 
 - A mobile device hotspot will work if WiFi is not available
 - Once a unit is activated, WiFi is not required for use
- Step by step instructions to activate the Mobilize Rescue App on Comprehensive & Mobile units can be found [HERE](#)
- Install tamper-proof seals and inspection card
Instructions: [COMPREHENSIVE](#) | [MOBILE](#)
- Additional materials can be found on the Mobilize Support Page [HERE](#)



PUBLIC ACCESS, UTILITY & COMPACT UNITS:

- Distribute customized Licensing Code to everyone in your community
- Instructions to activate the Mobilize Rescue App on smart phones can be found [HERE](#)
 - “Allow Notifications” for the Mobilize Rescue App to work properly
 - Please contact mobilize@zoll.com to request a customized licensing code
- Additional materials can be found on the Mobilize Support Page [HERE](#)



3. How to maintain the Mobilize Rescue Unit



Comprehensive & Mobile units need to be inspected and charged once a month and after any use

- Charge your unit by plugging it into a USB wall cube
- OR -
- Charge using the included portable charger
 - We recommend leaving the portable charger “full” between charges as a secondary power source
- Inspect the unit by following the Monthly Diagnostic Instructions: [COMPREHENSIVE](#) | [MOBILE](#)



Public Access, Utility & Compact units need to be inspected regularly and after any use

- Check that all medical supplies are intact, not expired and free of contamination
- Check for expired supplies

4. How to update the Mobilize Rescue App



- You will receive an email notification when updates are available
- WiFi or a personal hotspot will be needed for updates
- Update instructions for Comprehensive and Mobile units can be found [HERE](#)



- You will receive an email when smartphone app updates are available
- The Mobilize Rescue App updates are completed via your smart phones normal app update process

Why is it important to update the Mobilize Rescue App periodically?

- To address changes in medical guidelines
- To address changes in technology and/or operating systems

Go to www.mobilizerescue.com/support for more detailed instructions.
Questions? Email mobilize@zoll.com