

ZOLL



AED


ZOLL AED 3™

The AED that leads the way

SMART

Real CPR Help®

- Real CPR Help guides rescuers in delivering high-quality CPR consistent with current guidelines
- Real-time feedback on depth and rate of compressions
- Full-color images and integrated voice, text, and visual prompts for support throughout the rescue
- Enhanced bar gauge allows rescuers to see when they are doing good CPR in real-time



READY

Integrated Pediatric Rescue

- Universal CPR Uni-padz™ are always ready for both adult and child rescues

Low Cost of Ownership

- Electrode pads and batteries last five years – fewer replacements and a lower total cost of ownership

WiFi Connectivity

- Can automatically report AED readiness
- Quickly transfers case data

To learn more about the ZOLL AED 3™, call us at 1-(800) 804-4356 or visit us at www.zoll.com/aed3.

ZOLL®