

The Powerheart® AED G5 Trainer

Our simple-to-use trainer offers a realistic, hands-on rescue experience.

The Powerheart AED G5 Trainer plays a key role in your successful AED program.

Being rescue ready is more than just having AEDs available, it requires preparation of potential rescuers so they know what to do in a sudden cardiac arrest emergency. The Powerheart G5 Trainer was made to do exactly that.

By simulating different situations, students get valuable hands-on experience to handle a variety of rescue scenarios.

- Instructors can vary rescue conditions to test student responses using the wireless remote
- Scenarios can be paused so instruction can be done throughout the training
- The trainer will mimic automatic and semi-automatic AED operation
- Students can practice with reusable adult pads (with or without ICPR) and/or reusable pediatric pads

The Trainer provides an easy way to stay in practice.

Owning a G5 Trainer allows you to keep skills fresh and up to date to local regulations. It also gives you the flexibility to train as needed so that new employees or colleagues don't have to wait for an outside class to be rescue ready.

Users

- AED/CPR instructors
- Designated workplace responders
- EMS responders
- Anyone wanting to learn how to use an AED

Primary Benefits

Practice. Students gain real world experience to help a sudden cardiac arrest victim. And with CPR feedback, they get a good feel for effective CPR.

Ease of use. Clear, instructive voice prompts guide students through a simulated rescue, and the built-in metronome sets the proper chest compression pace. Configurable for use in a variety of languages for added flexibility.

Preparation. With four pre-programmed rescue scenarios, students will learn what to do in different situations.



The Powerheart® AED G5 Trainer



Overview	
Interactive	<ul style="list-style-type: none">· Four pre-programmed rescue scenarios following American Heart Association (AHA) and European Resuscitation Council (ERC) guidelines. Also follows the AHA Heartsaver AED curriculum.· Advance rescue stages and change parameters using an infrared, wireless remote.· Pause scenarios for continuous instruction throughout the training.· Voice prompts accurately represent what is heard during an actual rescue situation.· CPR feedback provides direction so students meet AHA/ERC CPR guidelines
Versatile	<ul style="list-style-type: none">· Change rescue scenarios to introduce different variables· Simulate various cardiac rhythms, both shockable and non-shockable· Monitor student skills during simulations· Configure the trainer with different languages as needed
Economical	<ul style="list-style-type: none">· Reusable training electrodes.· Compatible with any type of CPR manikin.· Trainer is powered with three C-cell alkaline batteries; wireless remote uses two AAA alkaline batteries.· Dual modes can be easily configured for semi-automatic and fully-automatic training.
Safety	<ul style="list-style-type: none">· Simulated delivery of defibrillation pulse for training purposes.· Cardiac Science AED Trainers will not generate or deliver defibrillation shocks.
Dimensions	<ul style="list-style-type: none">· 8 cm x 32 cm x 24 cm (3.3 in x 12.5 in x 9.5 in)
Part Numbers	<p>190-5020-001: G5 Trainer, Adult Electrodes 190-5020-002: G5 Trainer, Adult Electrodes with ICPR XTRAPAD004A: Reusable Adult Training Pads XTRAPAD005A: Reusable Pediatric Training Pads XTRAPAD006A: Reusable Adult ICPR Training Pad</p> <p>Package includes: Powerheart G5 Trainer, (1) Set of Adult Trainer Pads, (1) Infra-red wireless remote, (1) User's Guide CD with Language Packs, Overlay for automatic version, and Language Configuration Instructions</p>

